



SPICA

Your wheels in rehab

- Empowering mobility and independence during recovery.
- Designed for a perfect fit post hip or spinal surgery.
- Range of sizes and customisation ensuring optimal comfort and support.



Discover More



Enhanced comfort and mobility post surgery

Manual attendant-controlled Seating and Mobility rental service: Offering flexible solutions for short-term mobility and caregiving needs.

Specialist support for children/young adults in hip cast, spinal cast or radial extension braces. Highly configurable; accommodates awkward positions, including those created by large post-operative plaster cast.

Client profile

- Suitable for individuals from 3 years old up to 6'2" and 110 kg (17 stone).
- Migrating hip joint.
- Loss of range of movement affecting posture and function.
- Flexor contractions affecting hip & knee flexion, posture and function.
- Trauma to lower limb resulting in bone fracture (leg, foot).
- Continuity of care: vital contribution to posture maintenance and pressure management for existing wheelchair users.



BOOK A FREE DEMO

Feature


- Suitable for indoor and outdoor use – enhancing mobility.
- 45-Degree Tilt in Space; customizable positioning for optimal comfort and support.
- Individually adjustable and removable leg splints; accommodate abduction and knee flexion.
- NEW Standard options at no extra cost; swing away arm rests, activity tray, pommel, footbox.
- Minimum rental period; 3 Weeks – also ideal for holidays or traveling abroad.
- Crash test certified for use in Wheelchair Accessible Vehicles (WAV); compliant with ISO 7176-19:2022 standards.

Benefits

- Facilitates early discharge; transition safely from hospital to the comfort and security of home.
- Unchanged; friendly and functional for parents and carers short term use.
- Mobility with full access; suitable for all environments including home, school, and on-the-go.
- Adjustable seating positions; Ideal for feeding, activities, play, personal care and relaxation.
- Adaptable, versatile, and easy to manoeuvre; designed for effortless navigation, even in tight spaces.
- Suitable for transportation; enables continued participation in social, educational, and therapy environments.

Chair Dimensions and Sizing Guide

SPECIFICATION	EXTRA SMALL	SMALL	MEDIUM	LARGE
Max User Mass	110kg	110kg	110kg	110kg
Seat Cushion Depth Range	30-35cm	35-40cm	40-45cm	45-50cm
Seat Width	30cm	35cm	40cm	45cm
Back Rest Height	40cm	45cm	51cm	57cm
Total Backrest Height Growth	5cm	5cm	5cm	5cm
Backrest Width	25cm	30cm	31.5cm	35cm
Leg Rest Angle Range	80°-170°	80°-170°	80°-170°	80°-170°
Leg Rest Length Range	21-48cm	21-48cm	21-48cm	21-48cm
Back Rest Recline Angle Range	90-180°	90-180°	90-180°	90-180°
Tilt in Space	+3° to -45°	+3° to -45°	+3° to -45°	+3° to -45°



XS		
A	36-44cm	(14-17 in)
B	21-43cm	(8-17 in)
C	20-30cm	(8-12 in)
D	15-26cm	(6-10 in)

S		
A	42-52cm	(17-20in)
B	21-43cm	(8-17 in)
C	25-35cm	(10-14 in)
D	20-31cm	(8-12in)

M		
A	48-60cm	(19-24 in)
B	21-43cm	(8-17 n)
C	30-40cm	(12-16 in)
D	25-36cm	(10-14 in)

L		
A	54-66cm	(21-26 in)
B	21-48cm	(8-19 in)
C	35-45cm	(14-18 in)
D	30-41cm	(12-16 in)

Wheelbase Width			
E	Standard	62cm	(25 in)
	Narrow	58cm	(23 in)



#thechuncway

Custom-made indoor and outdoor wheelchairs with built-in postural support for children and adults with moderate to complex needs.

Designed, manufactured and hand-assembled in Herefordshire, England (in 10 working days*).

Supporting your everyday needs – and passionate about it too.

*97% all standard chairs made in 2023



Discover More

Web: chunc.com | **Tel:** 01432 377 512 | **Email:** info@chunc.co.uk